



Clinical Validation of CTC Subtype Frequency to Prognosis Overall Survival (OS) in Metastatic Castrate Resistant Prostate Cancer (mCRPC) Patients

Ryan Dittamore¹, Jessica Louw¹, Nicole Schreiber², Ryon Graf¹, Adam Jendrisak¹, Ann Johnson¹, Lyndsey Dugan¹, Brigit McLaughlin², Glenn Heller², Martin Fleisher², Yipeng Wang¹, Dena Marrinucci¹, Howard I. Scher^{2,3}



EPIC SCIENCES
www.epicsciences.com

¹ Epic Sciences, Inc., San Diego, CA

² Sidney Kimmel Center for Prostate and Urologic Cancers, Memorial Sloan-Kettering Cancer Center, New York, NY


³ Department of Medicine, Weill Cornell Medical College, New York, NY

Background

- CellSearch® clinical sensitivity is limited in mCRPC by exhibiting low CTC counts in many patients despite poor outcomes.
- CellSearch® detects a narrow phenotype of CTCs: EpCAM(+), CK(+), DAPI(+), CD45(+); which could reduce detection sensitivities.
- The Epic Sciences platform does not use enrichment, and detects an expanded range of CTC histology, including CK(+) CTCs, CK(-) CTCs, CTC Clusters, and Apoptotic CTCs.
- Expanded CTC subtypes require clinical validation as individual features.
- To assess the clinical value of CTC subtypes, CTC burden of expanded CTC subtypes was associated to overall survival and evaluated in a large cross-sectional cohort of modern mCRPC practice.

Methods

- 221 blood samples from 179 unique patients were collected prior to initiating Androgen Receptor (AR) directed (n = 150) or taxane (n = 71) therapy for mCRPC.
- Samples were analyzed with the Epic Sciences platform to enumerate CK(+) CTCs, CK(-) CTCs, CTC Clusters, Apoptotic CTCs, and Small CTCs (Figure A).
- Patients were followed for up to 2.3 yrs.
- Paired CellSearch® blood draws were processed at MSKCC Clinical Laboratory per manufacturer recommendations. CellSearch® counts were capped at 200 CTCs per tube (from 7.5 mL of blood). For comparison, CellSearch® and Epic Sciences counts were normalized per milliliter, capped at 26.7/mL. Paired CellSearch® and Epic Sciences traditional CTC counts were collected from 173 patient samples.



Schematic of Epic Sciences CTC Platform CTC enumeration, morphology, biomarker, & FISH analyses workflow:


- Nucleated cells from blood sample placed onto slides and stored in a -80°C biorepository
- Slides stained with cytokeratin (CK), CD45, DAPI, AR N-Term
- Slides scanned
- CTC candidates detected by a multi-parametric digital pathology algorithm
- Human reader confirmation of CTCs & quantitation of biomarker expression

Patient Demographics

Patient Primary Therapy	
Number of Unique Patients	179
Age, years	68 (45 – 91)
Primary Treatment	
Prostatectomy	84 (47%)
Radiation	34 (19%)
Brachytherapy	7 (4%)
None	54 (30%)
Characteristic All Samples	
Number of Baseline Samples	221
Age, years	68 (45 – 91)
Prior Hormone Therapies*	
1 - 2 lines	81 (37%)
3 lines	46 (21%)
≥4 lines	94 (42%)
Chemotherapy Status	
Chemo-naïve	136 (62%)
Chemo-exposed	85 (38%)
Metastatic Disease	
Bone Only	63 (29%)
Lymph. Node (LN) Only**	24 (10%)
Bone & LN	77 (35%)
Bone & Visceral + LN**	35 (16%)
Laboratory Measures	
PSA, ng/mL	37.7 (0.1 – 3728.2)
Hgb, (g/dl)	12.0 (7.0 – 15.0)
ALK, (unit/L)	110 (25 – 2170)
LDH, (unit/L)	222.5 (123 – 1293)
ALB, (g/dl)	4.2 (3.1 – 4.9)


Patient Line of Therapy			
Total Samples	1 st Treatment Decision	2 nd Treatment Decision	3 rd + Treatment Decision
A or E Baseline Blood Draw (n=150)	No Prior A or E (1st Line) n=64	Previous A or E (2nd Line) n=36	Previous A & E (3 rd + Line) n=5
T Baseline Blood Draw (n=71)	No Prior A or E (1st Line) n=12	Previous A or E (2nd Line) n=12	Previous A & E (3 rd + Line) n=10 Previous AR Tx & T (3 rd + Line) n=37

CTC Histological Subtypes Detected on the Epic Sciences CTC Platform




- CTCs enumerated in this study encompass several histological types.
- “Traditional” CK(+) CTCs are detected as single cells positive for cytokeratin expression (Figures B-C).
 - Some CK(+) CTCs are smaller than surrounding white blood cells (Figure C).
 - CK(-) CTCs have distinctive nuclear malignant features and/or the presence of AR N-terminal domain (Figures D-E).
 - CTC Clusters consist of more than one adjacent CTC (Figures F-G).
 - Apoptotic CTCs (Figures H-I) contain fragmented nuclei.
 - AR localization can be cytoplasmic (Figure B), nuclear (Figures D-F, H-I) or both, even within a single CTC Cluster (Figure G).

Clinical Sensitivity of CTC Detection: CellSearch® vs. Epic Sciences



- 173 of 221 patient samples had matched CellSearch® testing.
- “Epic” is the sum of all CTC subtypes detected on the Epic Sciences platform per sample: CK(+), CK(-), CTC Clusters, Apoptotic, and Small CTCs.
- CellSearch® vs. Epic Sciences enumeration in matched samples shown side-by-side in a matched bar plot (Figure J).
- The difference between Epic Sciences and CellSearch® enumeration is shown per sample (Figure K).
- Note: CellSearch® counts were capped at 200 per tube (from 7.5 mL of blood) by MSKCC clinical laboratory. For comparison, CellSearch® and Epic Sciences counts were normalized per milliliter, capped at 26.7/mL



Platform Comparison of CTC Detection


(L)	All Patients (n = 173)	1st line (n = 55)	2nd line (n = 44)	3rd line (n = 34)	4 th + line (n = 40)
CellSearch®: CTCs Present	116 (67%)	32 (58%)	24 (55%)	27 (79%)	33 (83%)
CellSearch®: CTCs Absent	57 (33%)	23 (42%)	20 (45%)	7 (21%)	7 (17%)
Epic Sciences: CTCs Present	166 (96%)	52 (95%)	43 (98%)	31 (91%)	40 (100%)
Epic Sciences: CTCs Absent	7 (4%)	3 (5%)	1 (2%)	3 (9%)	0 (0%)

(M)	All Matched Samples (n = 173)	(N)	CellSearch® Favorable Count (< 5 CTC / 7.5mL = <0.67/1mL, n = 102)	
	CellSearch®	Epic Sciences	CellSearch®	Epic Sciences
Median	0.27 / mL	6.82 / mL	0 / mL	4.37 / mL
Range	0 - 26.7* / mL	0 - 991 / mL	0 - 1 / mL	0 - 144 / mL

- CTC detection by line of therapy in mCRPC (Figure L).
- Median and range of CellSearch® and Epic CTC detection (Figure M).
- Enumeration in samples with CellSearch® count less than prognostic threshold of 5 CTCs/tube (0.67/mL) (Figure N).
- Note: Epic tested two slides (corresponding to roughly 1mL of blood) while CellSearch® assayed 7.5mL, potential bias against Epic in comparison.

*CellSearch® and Epic counts were normalized per milliliter, capped at 26.7/mL

CTC Histological Subtypes Prognosticate Overall Survival



Univariate thresholds for KM plots were chosen based on time-dependent survival ROC curves. CTC/mL was included in multivariate models as a log2-transformed continuous variable as previously reported in Scher et al 2009.

Conclusions

- The Epic Sciences platform has increased clinical sensitivity for mCRPC CTC detection rate vs. CellSearch®, (96% vs. 67%) and magnitude of enumeration (median 6.82/mL vs. 0.27/mL)
- All subtypes of CTCs detected by Epic Sciences: CK(+), CK(-), Small, Apoptotic, Clusters are prognosticators of shorter OS in univariate models.
- CTC Subtypes: CK(+), CK(-) and Small CTCs as well as all Epic Sciences CTC Subtypes pooled, each add to the prognostication of OS in multivariate models
- Characterization of non-traditional CTCs (CK- and Small CTCs) provides increased clinical sensitivity and may provide key insights to cancer biology.

Support: MSKCC SPORE in Prostate Cancer (P50 CA92629), the Department of Defense Prostate Cancer Research Program (PC051382), The Prostate Cancer Foundation, Mr. William H. Goodwin and Mrs. Alice Goodwin and the Commonwealth Foundation for Cancer Research, The Experimental Therapeutics Center of Memorial Sloan-Kettering Cancer Center.