


Background


The NCI 9012 study, which randomized 1st line mCRPC patients to Abiraterone or Abiraterone + Veliparib, demonstrated no advantage for Veliparib in any pre-specified subgroup, with DNA damage response mutations (DDRM) pts having better outcomes on both arms. This contrasted with the TOPARP trial¹, which enrolled heavily pre-treated pts with high CTC counts, where responses to PARPi were observed, particularly in patients with DDRm. Here we explore the use of CTCs in the NCI 9012 study to ascertain phenotypic Genomic Instability (pGI) as a proposed selection biomarker for PARPi.

The pGI CTC Biomarker Methodology


I - Schematic of the pGI (phenotypic Genomic Instability) algorithm training and validation



II - Example CTC images and corresponding single cell sequencing from the training set



III - Prognostic value of pGI status at Baseline (mCRPC ARSI treated pts)



- The pGI algorithm takes as an input the individual CTC image and predicts whether or not the CTC has a high number of Large Scale Transitions (i.e., Genomically Instable). The algorithm was trained on the ground truth determined by single cell Next Generation Sequencing (NGS).
- Example CTC images and corresponding single cell copy number variation (CNV) plots by chromosome for 3 representative CTCs from the training set. CTC-1 has a 'flat' or unperturbed genome while CTC-2 and CTC-3 have a high number of LSTs and are therefore, genomically unstable.
- Independent testing of pGI in a cohort of mCRPC Pts treated with ARSI Rx (Abiraterone or Enzalutamide), demonstrating that pGI is a strong prognostic indicator of poor OS.

Patient demographics and study design

Patient Characteristic	Arm A (Abi) N=35	Arm B (Abi + Veliparib) N=40	Significance	p-value*
Median Age in yr (min, max)	67.2 (50.6, 79.8)	68.1 (53.2, 85.8)	NS	0.85
CTCs > 3	29 (82.86%)	32 (80.00%)	NS	1.00
Median CTC/mL (min, max)	2.8 (0, 259)	3.8 (0, 287)	NS	0.55
pGI+ CTC/mL > 0	12 (34.29%)	23 (57.50%)	NS	0.07
Median pGI+ CTC/mL (min, max)	0 (0, 129)	1.0 (0, 199)	NS	0.37
Race - no. (%)			NS	0.50
White	32 (91.43%)	37 (92.50%)		
Black	2 (5.71%)	1 (2.50%)		
Other	1 (2.86%)	2 (5.00%)		
ECOG - no. (%)			NS	0.32
0	22 (62.86%)	29 (72.50%)		
1	13 (37.14%)	10 (25.00%)		
2	0 (0.00%)	1 (2.50%)		
Median PSA ng/mL (min, max)	28.94 (0.8, 458.0)	37.32 (0.04, 802.0)	NS	0.24
Cancer Pain Present - no. (%)	12 (34.29%)	10 (25.00%)	NS	0.58
Sites of Disease - no. (%)			NS	0.25
Lymph Nodes	21 (60.00%)	28 (70.00%)		
Bone	32 (91.43%)	35 (87.50%)		
Visceral	8 (22.86%)	11 (27.50%)		
Other	6 (17.14%)	9 (22.50%)		
Previous Treatments - no. (%)			NS	0.21
Chemotherapy	5 (14.29%)	12 (30.00%)		
Enzalutamide	1 (2.86%)	0 (0.00%)		
Experimental Agent	9 (25.71%)	4 (10.00%)		
Surgery T	30 (85.71%)	7 (17.50%)		
ETS Fusion - no. (%)			NS	1.00
Positive	14 (40.00%)	16 (40.00%)		
Negative	21 (60.00%)	24 (60.00%)		
Previous Metastases - no. (%)	3 (8.57%)	4 (10.00%)	NS	1.00
Median Treatment Cycles (min, max)	9 (3, 39)	12 (3, 38)	NS	0.46
Confirmed PSA Response - no. (%)	22 (62.86%)	30 (75.00%)	NS	0.32
Confirmed PSA Progression - no. (%)	25 (71.43%)	29 (72.50%)	NS	0.84
Median Time to PSA Progression (Months)	8.57 (0.89, 25.3)	11.03 (0.89, 25.3)	NS	0.27
Measurable disease at baseline - no. (%)	20 (57.14%)	25 (62.50%)	NS	0.49
Median PSA (Months)	8.28 (0.89, 25.3)	11.03 (0.89, 25.3)	NS	0.44

*p-values calculated using Fisher's exact test for categorical variables and Welch's t-test for continuous variables

82 Pts from NCI 9012
1717 CTCs identified at Baseline (81.3% of Pts)
1073 CTCs predicted GI+ at Baseline (46.6% of patients)
31 Matched tumor biopsy sequencing for DNA repair deficiencies at baseline


5 Pts Clinically Unevaluable
2 Pts Failed QC

Arm A: Abiraterone
29/35 (82.86%) Pts CTC+
12/35 (34.29%) Pts pGI+
5/13 (38.46%) Tissue DDRm+


Arm B: Abiraterone + Veliparib
32/40 (80.00%) Pts CTC+
23/40 (57.50%) Pts pGI+
5/18 (27.78%) Tissue DDRm+

Endpoint: PSA Progression

pGI+ and total CTC prevalence




pGI positive patients associate with improved response in the Abi + Veliparib Tx arm vs. Abi alone




A) PSA waterfall. B) Kaplan-Meier analysis and Fisher's exact test comparing PSA responders and non-responders by Tx arm. Any patient having greater than 0 pGI+ CTCs is defined as pGI positive.


CTC pGI and Tissue DDRm relationship to outcome




pGI+ CTCs decrease after start of Abi + Veliparib but not Abi alone




Patient Examples of Single CTC Genomic Instability and Clonal Heterogeneity & Relationship to Outcome



DDR Genomic alterations observed in CTC before and at progression of Abi + Veliparib Tx




Genomic Instability & DDR alterations observed in CTCs in patient with non-response to Abi



Discussion

- Phenotypic Genomic Instability (pGI) CTC algorithm identifies single cell and patients with true genomic instability
- Patients with pGI+ CTCs receiving Abiraterone + Veliparib vs. Abiraterone alone are associated with improved PSA response
- pGI CTCs have a significant reduction from baseline to Post-Tx-start in Abiraterone + Veliparib but not in Abiraterone only
- Single CTC genomic analysis supports biologic mechanism of response & resistance to Veliparib & Abiraterone

Analytical performance of the pGI biomarker



References
1. Kato, Joaquin, et al. "DNA-repair defects and olaparib in metastatic prostate cancer." *N Engl J Med* 373.18 (2015): 1637-1708.
2. Scher, Howard I., et al. "CTC phenotype classifier to identify mCRPC patients (pts) with high genomic instability CTCs and to predict failure of androgen receptor signaling (AR Tx) and taxane (T) systemic therapies." (2016): 5044-5044.