



# Integrated Single Cell Phenogenomic Subtyping of CTCs Identify Inter-Cellular Tumor Heterogeneity (het) and Multiple Resistance Mechanisms in mCRPC Patients (pts)

Howard I. Scher<sup>1</sup>, Ryon Graf<sup>2</sup>, Adam Jendrisak<sup>2</sup>, Nicole A. Schreiber<sup>1</sup>, Brigit McLaughlin<sup>1</sup>, Stephanie Green<sup>2</sup>, Angel Rodriguez<sup>2</sup>, Martin Fleisher<sup>1</sup>, Jerry Lee<sup>2</sup>, James Kelvin<sup>2</sup>, Mahipal Suraneni<sup>2</sup>, Yipeng Wang<sup>2</sup>, Mark Landers<sup>2</sup>, Ryan Dittamore<sup>2</sup>

EPIC SCIENCES  
www.epicsciences.com

<sup>1</sup> Sidney Kimmel Center for Prostate and Urologic Cancers, Memorial Sloan-Kettering Cancer Center, New York, NY  
<sup>2</sup> Epic Sciences, Inc., San Diego, CA




## Background


Intra-tumor and inter-cellular heterogeneity is highly prevalent in progressing mCRPC pts and has been associated with resistance to ARSi Tx.<sup>1</sup> While genomic profiling of metastatic tumors may identify pts with actionable targets, the response to novel Tx outside of ARSi, taxanes and DDR directed drugs are infrequent and not durable. Basket trials have also had limited impact for mCRPC pts. We sought to better understand tumor heterogeneity at the single cell level for the purposes of:

- Better defining the sensitivity of the individual circulating tumor cell (CTC) subtypes to ARSi Tx & taxane-based chemotherapy
- Better understanding tumor evolution and better characterization of subclonal tumor types, mechanisms of resistance & clustering of genetic changes
- Developing more effective biomarker directed treatment approaches for heterogeneous tumors


## CTC Based Phenogenomic Analysis Using Digital Pathology and NGS Reveals CTC SubTypes Associated with Clinical Outcomes




## Patient Demographics




## Presence of Phenotypic CTC Subtypes Prognosticate Worse Outcomes on ARSi & Taxane Therapies



## PhenoGenomic Profiles of CTC Subtypes




## Phenotypic CTC Subtypes are Associated with Distinct Gene Signatures



A. Representative cell images of 15 phenotypic CTC subtypes.  
B. Distribution of large scale transitions (LSTs), a surrogate of genomic instability. LST was measured as *n* of chromosomal breaks between adjacent regions of at least 10 Mb.  
C. Prevalence of 5 representative cancer gene copy number alterations.

## Specific CTC Subtypes are Enriched in Samples with High Heterogeneity Scores



## Conclusions

- Detecting specific CTC phenotypic subtypes in patients predict the relative risk of failure to ARSi and taxanes, and may inform Rx selection
- CTC subtypes (E, G, K, I, & M) are detected more frequently in high CTC heterogeneity samples (adjusted p-value = 0.05)
- Phenogenomic analysis of CTC subtypes that predict for resistance to ARSi and taxanes and high heterogeneity harbor novel genomic patterns associated with differential prognosis of survival. e.g.: The presence of Cell Type K, enriched for PTEN and p53 loss, is associated with worse prognosis on ARSi and the presence of Cell type G, enriched for RB1 loss showed worse prognosis on Taxanes
- Identifying different CTC subtypes may provide insight into disease evolution, inform treatment decisions for individual patients and guide new drug development. Longitudinal monitoring of phenotypic CTC subtypes in multiple clinical trials is currently underway

**Support:** MSKCC SPORE in Prostate Cancer (P50 CA92629), the Department of Defense Prostate Cancer Research Program (PC051382), The Prostate Cancer Foundation, Mr. William H. Goodwin and Mrs. Alice Goodwin and the Commonwealth Foundation for Cancer Research, The Experimental Therapeutics Center of Memorial Sloan-Kettering Cancer Center.