

Improved sensitivity in identification of ER and HER2 expressing metastatic breast cancers with a combination of cell & cell free liquid biopsy analysis

Giuseppe Di Caro¹, Ernest Lam¹, Kandra Horne¹, Stefan Gluck¹, Megan Slade¹, Rick Wenstrup¹, Lee Schwartzberg²

¹Epic Sciences, San Diego, CA, ²Medical Oncology and Hematology, Renown Institute for Cancer, University of Nevada, Reno, Nevada



BACKGROUND


Liquid biopsies are a non-invasive diagnostic approach for detecting circulating tumor cells (CTCs) or circulating tumor DNA (ctDNA) that may provide clinically actionable information to guide treatment decisions for metastatic breast cancer (MBC) patients when a conventional biopsy is otherwise infeasible. In addition, the development of HER2 and ER assays is expected to enable the identification of HER2 expressing MBC that may benefit from novel HER2-targeted therapies. Here we report a comprehensive liquid biopsy platform including immunofluorescent HER2 and ER protein expression in CTCs (ctCIF) coupled with the determination of ERBB2 amplification by single-cell CTC genomics (ctcDNA), and ctDNA alterations in plasma

METHODS

Blood samples from 247 progressing metastatic breast cancer patients were collected for cell & cell-free DNA analysis. Blood was also collected from 25 blood donors (HDs) with no known cancer history as controls. After plasma isolation, nucleated cells were plated, & slides were bio-banked. Immunofluorescent staining & subsequent imaging were performed on replicate slides. CTCs were identified using Epic Sciences digital imaging & machine learning algorithms. Single-cell CTC isolation and sequencing for genomic quantification of large-scale state transitions (LSTs) and ERBB2 Copy Number Variants (CNV) was used to determine genomic instability and ERBB2^{amp} in individual CTCs. Bio-banked plasma was analyzed to detect ctDNA alterations with high clinical relevance (Class IA).


EPIC SCIENCES LIQUID BIOPSY PLATFORM

Figure 1. CTC morphology by digital pathology coupled with Immunofluorescent detection & biomarkers expression (HER2⁺ and ER⁺ CTCIF)



Schematic of Epic Science's CTC collection and detection process: 1) Nucleated cell from blood are plated onto slides and stored in the -80C biorepository; 2) Slides are stained with CK, CD45, DAPI and HER2; 3) Slides are scanned, and CTC candidates detected by multi-parametric digital pathology algorithm.


Figure 2. CTC single cell isolation for low-pass Whole Genome Sequencing (WGS) (LST⁺ and ERBB2^{amp} CTCDNA)



Schematic of Epic Science's CTC single cell genomics: 1) CTC are picked with a microcapillary and isolated for single cells WGS 2) Large scale transitions (LST) values and ERBB2^{amp} Z scores are computed through the bioinformatics pipeline for each CTC candidate.

Figure 3. Plasma isolation & ctDNA analysis for SNVs and CNVs detection

56 gene ctDNA panel, TMB, MSI



Epic Science's ctDNA pipeline: Plasma is isolated and SNVs Variant Allele Frequencies (VAF) and genomics CNVs relevant for breast cancer are computed through the bioinformatics pipeline for each patient.

RESULTS (CTC⁺DNA AND CTC⁻DNA GENOMIC, AND CTC⁺IF PLATFORMS)

Prevalence of CTCs, CTC with genomic instability (CTC LST⁺), ctDNA (all, and 1A) and ERBB2^{amp} by ctDNA and CTCDNA genomics

	CTC ⁺ %	CTCDNA (LST ⁺) %	ctDNA (all*) %#	ctDNA (IA*) %#	ERBB2 ^{amp}	
					ctDNA(IA*) %#*	CTCDNA %#^
MBC (n=247)	94 (233/247)	63 (157/247)	62 (112/182)	30 (54/182)	2 (3/182)	15 (37/247)
Healthy donors (n=25)	0 (0/25)	0 (0/25)	0 (0/25)	0 (0/25)	0 (0/25)	0 (0/25)

65 patients had failed QC for ctDNA analysis (or NSQ non-sufficient ctDNA quantity)

*ctDNA LoD ≥ 1.4 "apparent" folds change (VAF/CNV) ^CTCDNA LoD ≥ 2 folds change (single cell level)

Table 1 (Left). No false positive detection in Healthy Donors, suggesting high specificity for each analysis. Table 1 (Right). CTCDNA single-cell genomics was more sensitive than ctDNA in detecting ERBB2^{amp} in MBC patients (15%, and 2%, respectively).

Distributions of HER2 and ER protein MFI in WBC and CTC by CTCIF




Figure 4. HER2 (left) and ER (right) MFI expression by CTCIF across all patients slides in White Blood Cells (WBC), and CTC identify HER2 (left) and ER (right) negative, and low threshold values.

Prevalence of HER2⁺ and ER⁺ CTC

	CTCIF HER2 ^{exp} %	CTCIF ER ^{exp} %
MBC (n=247)	79 (194/247)	76 (189/247)

Table 2. HER2⁺ and ER⁺ CTC by CTCIF were detected in 79% and 76% of MBC patients, respectively.

Liquid biopsy classification of HER2 status

	HER2 amplified <i>ERBB2</i> ^{amp} by CTCDNA or ctDNA %	HER2 expressed HER2 non-amplified HER2 ^{exp} by CTCIF and <i>ERBB2</i> ^{non-amp} by CTCDNA or ctDNA %	HER2 negative HER2 ^{non-exp} by CTCIF and <i>ERBB2</i> ^{non-amp} by CTCDNA or ctDNA %
MBC (n=247)	15 (37/247)	65 (162/247)	20 (48/247)

Table 3. A liquid biopsy classification of HER2 status by combining the three platforms (CTCIF, CTCDNA, and ctDNA) identified that among MBC, 15% were ERBB2^{amp}, 65% were HER2^{exp} (HER2⁺ and ERBB2^{nonamp}), and 20% were HER2^{neg} (HER2^{nonexp} and ERBB2^{nonamp}).

Liquid biopsy classification of combined HER2 and ER status in ERBB2 non amplified MBC




Table 4. A liquid biopsy classification of combined HER2 and ER status by combining the three platforms (CTCIF, CTCDNA, and ctDNA) identified that among ERBB2^{nonamp} MBC: 58% were ER^{exp}, and HER2^{exp}; 15% were ER^{nonexp}, and HER2^{exp}; 15% were ER^{exp}, and HER2^{nonexp}; 8% were ER^{nonexp}, and HER2^{nonexp}.

CONCLUSIONS

1. This liquid biopsy platform integrates comprehensive CTC analysis and ctDNA to report HER2 and ER status with high specificity and sensitivity in MBC patients.
2. CTCDNA single-cell genomics is more sensitive than ctDNA at detecting ERBB2^{amp}.
3. This platform has potential clinical utility in selecting appropriate treatments for MBC without a tissue biopsy and/or in later lines of therapy among often biologically heterogeneous tumor sites of metastatic disease.