- conventional tissue biopsies are not feasible.
- metastatic triple-negative breast cancer (mTNBC), which has limited treatment options.
- patients, including patients with low levels of each biomarker.





expression levels which follow the same trend demonstrated with RNA-Seq data. Three of these cell lines (highlighted in the square boxes) were selected for use during development of the TROP2 assay.

TROP2 and HER2 Expression by Liquid Biopsy in Women with Metastatic Triple Negative Breast Cancer Alessandra Cunsolo¹, David Bourdon¹, Brandon Guillory¹, Megan Slade¹, Giuseppe Di Caro¹, Alisa Tubbs¹, Naoto T. Ueno² ¹Epic Sciences, San Diego, CA. ²Department of Breast Medical Oncology, The University of Texas MD Anderson Cancer Center, Houston, TX, USA

a single TROP2⁽⁺⁾ CTC. Bottom row, a cluster of two distinct TROP2⁽⁺⁾ CTCs.



TROP2 and HER2 IF Assay Signals in mTNBC Patient Samples



18 mTNBC Patient Samples Stained with Epic Sciences TROP2 and HER2 Assays



Figure 6: Assessment of TROP2 and HER2 expression among 18 mTNBC samples.

89% of patients (16 patients) with mTNBC had detectable CTCs. (A) Each dot on the plot is an individual CTC from a patient. TROP2 MFI, and HER2 MFI across all cells had a wide dynamic range (mean:1357, range: 37-38281, and mean:634, range: 35-24985, respectively). (B) TROP2 and HER2 positivity proportions within patients with CTCs.

CONCLUSIONS

- We report on developing a TROP2 immunofluorescent assay that utilizes Epic Sciences' CTC detection platform.
- Across exploratory mTNBC patients, a high prevalence (89%) of TROP2 expressing CTCs was detected with a wide dynamic range in biomarker signal.
- Blood-based assessment of mTNBC patients could further characterize relevant cancer phenotypes and may provide additional information for decision-making for anti-TROP2-antibody conjugate drugs (ADC) based on their mechanisms of action.

This presentation is the intellectual property of the author/presenter. Contact them at alessandra.cunsolo@epicsciences.com for permission to reprint and/or distribute.

